Red Blood Cell Lysis Buffer 1x

Ordering info

TBB0400, RBC Lysis Buffer 1x, 100 mL

Description

Red Blood Cell Lysis Buffer 1x is an optimized buffer designed for the lysis of human and murine red blood cells, with minimal effects on white blood cells. It allows to use leucocytes for further applications (cytometry, DNA/ RNA extraction).

Features

- Ready to use.
- Sterile.
- **Efficient and quick lysis** of red blood cells.
- **Cells remain viable** after red blood cell lysis.
- pH=7.1-7.4.

Applications

- Isolation of peripheral blood mononuclear cells.
- Isolation of T-cells for cytometry.
- Isolation of B-cells.
- Isolation of DNA, RNA.
- Immunofluorescence staining of whole blood for cytometry.

Storage

Store at 2-8 °C.

Quality Control

- pH tested.
- Lysis of human peripheral blood (A550 nm < 0.3).

Material required (not supplied)

- Water
- PBS 1x pH 7.4

Also available:

WATER, nuclease free (TBB0300, TBB0301)

PBS 1x pH 7.4 (TBB0360, TBB0361)

PROTOCOL

Warm RBC Lysis Buffer 1x to room temperature prior to use!.

Lysis of Human Peripheral Red Blood Cells

1. Add **2 mL RBC Lysis Buffer 1x** to each tube containing **up to 100 µL anti-coagulated whole blood**.
2. Immediately vortex the sample.
3. Incubate at room temperature, 10-15 minutes. Protect the reaction from light.
4. Centrifuge 350-500 g, 5 minutes. Discard the supernatant.
5. Wash the pellet with **1 mL of appropriate buffer** or **PBS 1x pH 7.4**.
6. Proceed with further procedures.

Lysis of Mouse Spleen Red Blood Cells

1. Cut with a scissors the harvested mouse spleen in small pieces.
2. On a petri dish, add the pieces in **5 mL complete media** or **PBS 1x pH=7.4**.
3. Mash the tissue with a plunger until an homogeneous mixture is obtained.
Add 5 mL complete media or PBS extra if is necessary.
4. Centrifugate at 350-500 g, 10 minutes and discard the supernatant.
5. Resuspend the pellet in **5 mL RBC Lysis Buffer 1x**.
6. Incubate on ice, 5 minutes with occasional shaking.
7. Add **30 mL PBS 1x pH 7.4**.
8. Centrifugate at 350-500 g, 5 minutes. Discard the supernatant.
9. Wash the pellet with **1 mL of appropriate buffer** or **PBS 1x pH 7.4**.
10. Proceed with further procedures.

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