Saliva Genomic DNA Purification Kit

Ordering info

TBK0145, 3 reactions (sample)

TBK0146, 50 reactions

TBK0147, 200 reactions

Description

Saliva Genomic DNA Purification Kit is an optimized kit based on salting-out principle to produce higher quantity and quality DNA. The kit is suitable for isolate genomic DNA from human and animal saliva.

Features

- High yield and purity, 1-10 µg, A₂₆₀/A₂₈₀ ~1.8.
- Scalable, easily to process many samples simultaneously.
- Safe, no phenol extraction.
- Fast, easy and cost-effective protocol.

Applications

DNA obtained is suitable for downstream molecular biology applications such as PCR, enzymatic digestion for cloning or Southern, genotyping, etc.

Quality Control

DNA isolation from saliva is checked by: integrity (agarose gel electrophoresis), quantity and quality (A₂₆₀/A₂₈₀ = 1.8 ± 0.2).

Kit Components

Components	TBK0146	TBK0147
PBS 1x pH=7.4	60 mL	250 mL
BS2 Buffer	12 mL	50 mL
BS3 Buffer	5 mL	20 mL
BS4 Buffer	0.6 mL	2 x 1.1 mL
Proteinase K	20 mg ^a	4 x 20 mg ^a
Proteinase K Resuspension Buffer	1 mL	4 x 1 mL
Elution Buffer	15 mL	25 mL

Order Info Kit Components: PBS 1x pH 7.4 (TBB0360) | BS2 Buffer (TBB0507) | BS3 Buffer (TBB0508) | BS4 Buffer (TBB0509) | Proteinase K (TBZ0305) | Proteinase K Resuspension Buffer (TBB0546) | Elution Buffer (TBB0510).

Before its use:

- ^aAdd 1 mL Proteinase K Resuspension Buffer and mix well.
- Store at -20°C.

Storage

- Store the kit at 25°C.
- Store BS4 Buffer at -20°C.
- Store Proteinase K at 2-8°C (short storage) or at -20°C (long storage).

Material required (not supplied)

- Sterile container.
- 2 mL Centrifuge Tube.
- Isopropyl Alcohol (CAS 67-63-0).
- Ethanol 70%.


PROTOCOL

I. SALIVA SAMPLE COLLECTION

- Transfer **1 mL Saliva** in 2 mL tube.
Spit several times into a sterile container. Avoid to drag nasopharyngeal secretions.
- Add **1 mL PBS pH 7.4** and mix well.
- Centrifuge at **13,000 g** for **3 minutes**. Remove the supernatant leaving a residual volume of **~50 µL**.
- Resuspend the cell pellet by vortex for **10-15 seconds**.

II. DNA PURIFICATION

- Add **20 µL Proteinase K** and **600 µL BS2 Buffer**. Mix by pipetting.
- Incubate at **55°C** for **0.5 - 1 hour**. Cool down to room temperature,
- Add **70 µL BS3 Buffer** and mix vigorously by vortex for **20 seconds**.
- Incubate at **4°C** for **5 minutes**.
- Centrifuge at **13,000 g** for **5 minutes**.
- Recover supernatant by transferring it to a clean 1.5 mL tube containing **190 µL Isopropanol** and **10 µL BS4 Buffer**.
- Mix by gentle inversion (**~ 50 times**).
- Centrifuge at **13,000 g** for **2 minutes** and discard carefully the supernatant using a pipette.
DNA is visible like a white pellet.
- Add **200 µL Ethanol 70%**.
- Centrifuge at **13,000 g** for **2 minutes** and remove the supernatant.
- Invert and dry the tubes on clean paper towels for **5-10 minutes**.
- Add **50-100 µL Elution Buffer** or **Water, Molecular Biology Grade** and resuspend DNA isolated.
- Check DNA quality on agarose electrophoresis gel and quantity by spectrophotometry. Store at **-20°C**.

Saliva Genomic DNA Purification Kit

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Componentes Kit: PBS 1x pH 7.4 (TBB0360) | BS2 Buffer (TBB0507) | BS3 Buffer (TBB0508) | BS4 Buffer (TBB0509) | Proteinase K (TBZ0305) | Proteinase K Resuspension Buffer (TBB0546) | Elution Buffer (TBB0510).

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