TIARIZOL™ REAGENT

Ordering info

TBR0100, TIARIZOL™ Reagent, 50 mL

TBR0101, TIARIZOL™ Reagent, 100 mL (2 x 50 mL)

Description

TIARIZOL™ REAGENT is a red monophasic ready-to-use formulation widely used and proven for RNA and DNA isolation. It allows a quick, economical, and efficient isolation of total RNA/ DNA. The reagent is employed in a sequential procedure: cationic lysis, followed by organic extraction and alcohol precipitation. After lysis and homogenization of the biological sample, three phases are separated. RNA and DNA are precipitated with isopropyl alcohol from colorless aqueous phase (upper) and interphase, respectively.

Features

- Ready to use formulation.
- Easily scalable nucleic acid isolation.
- Simultaneous RNA and DNA extraction.
- High yield and purity (A260/A280 ~2.0).
- RNA obtained can be directly used for any downstream applications: cDNA library, Northern blotting, nuclease protection assay, RT-PCR, in vitro translation, etc.

Applications

- Isolation of total RNA and DNA from various types of tissues and cells.
- RNA and DNA isolation from samples of different sources: human, animal, plant, yeast, bacterial, viral or cell culture.

Storage

Store at 2-8°C. Protect from light.

Reagents and Solutions Required (not supplied)

- Chloroform (CAS 67-66-3).
- Isopropyl Alcohol (CAS 67-63-0).
- Ethanol (CAS 64-17-5).
- Sodium Citrate, Dihydrate (CAS 6132-04-3).
- Sodium Citrate 0.1M in 10% Ethanol, pH=8.5.
- Ethanol 70% in Water Nuclease Free.
- Water, nuclease free (TBB0300).

Quality Control

Functionally tested.

Also available:

High-Q Spin Column TIARIZOL™ Plus RNA Purification Kit (TBK0244, TBK0245)

TIARIS™ Protect Tissue Buffer (TBB0430)

PROTOCOL FOR TOTAL RNA ISOLATION

I. HOMOGENIZATION

Blood: Add **3 volumes of TIARIZOL™ Reagent** to 300 µL of Blood or Buffy Coat and mix well by vortex.

Tissues: Use the following proportion:

50-100 mg tissue	1 mL TIARIZOL™ Reagent
1-10 mg tissue	0.8 mL TIARIZOL™ Reagent

OPTIONAL: Following homogenization, insoluble material could be removed by centrifugation at 12.000 x g for 10 minutes at 4°C. Transfer the cleared homogenate to a fresh tube.

Cells grown in suspension: Pellet cells at 200 x g for 5 minutes at room temperature. Resuspend cell pellet using the following proportion:

5.10 ⁶ cells	1 mL TIARIZOL™ Reagent
10 ² -10 ⁶ cells	0.8 mL TIARIZOL™ Reagent

Pass the lysate several times through a pipette tip.

At his stage, samples can be stored for at least one month at -80°C.

Cells grown on monolayer: Lyse cells directly in a culture dish or flask by adding 1 mL of TIARIZOL™ Reagent per 10 cm² growth area, pipette the cell lysate several times to ensure sufficient cell disruption.

II. PHASE SEPARATION

1. Incubate 5 minutes at room temperature.
2. Add **0.2 mL of chloroform** per 1 mL of TIARIZOL™ Reagent used.
3. Shake vigorously by vortex for 15 seconds.
4. Incubate 3 minutes at room temperature.
5. Centrifuge at 12,000 g for 15 minutes (or 2,600 g for 30 minutes) at 4°C.
6. The sample will separate into a red lower organic layer, a white interphase and a colorless upper aqueous layer that contains the RNA.

III. RNA PRECIPITATION

1. Transfer the aqueous phase very carefully, without disturbing the interphase to another tube.
2. Add **0.5 mL of cold isopropyl alcohol** per 1 mL TIARIZOL™ Reagent used.
3. Incubate for 10 minutes at room temperature.
4. Centrifuge at 12,000 g for 10 minutes (or 2,600 g for 30 minutes) at 4°C.

IV. RNA WASH

1. Remove the supernatant and wash the pellet once with 70% ethanol, adding at least **1 mL Ethanol 70%** per 1 mL of TIARIZOL™ Reagent used.
2. Centrifuge at 7,500 g for 5 minutes at 4°C and remove the supernatant.

V. RNA RESUSPENSION

1. Air-dry the pellet and dissolve in Water (nuclease free) by pipetting the solution up and down.
2. Incubate for 10 minutes at 60°C if necessary.
3. Store RNA at -80°C.

PROTOCOL FOR DNA ISOLATION

I. HOMOGENIZATION

II. PHASE SEPARATION

Employ the same procedures described in previous protocol. DNA is partitioned in the white interphase.

III. DNA PRECIPITATION

1. Carefully remove all upper colorless aqueous phase.
2. Change pipet tip, transfer the white interphase to a clean tube.
3. Add 0.3 mL of absolute ethanol per 1 mL of TIARIZOL™ Reagent used in the previous steps.
4. Mix well by inverting tubes several time.
5. Incubate for 10 minutes at room temperature.
6. Centrifuge at 2,000 g for 10 minutes at 4°C.

IV. DNA WASH

1. Remove the supernatant and add 1 mL of Sodium Citrate Solution per mL of TIARIZOL™ Reagent used at the beginning.
2. Incubate at room temperature for 30 minutes. Gently invert the tube. Repeat the inversion regularly.
3. Centrifuge at 2,000 g for 10 minutes at 4°C.
4. Repeat once steps 1-3.
5. Remove the supernatant and wash the pellet with 70% ethanol, adding 1.5 mL of ethanol per 1 mL of TIARIZOL™ Reagent used.
6. Centrifuge at 2,000 g for 5 minutes at 4°C.
7. Discard the supernatant.

V. DNA RESUSPENSION

1. Air-dry the pellet for 10 minutes and dissolve in Water (nuclease free) by pipetting.
2. Incubate for 10 minutes at 60°C if necessary.
3. Store DNA at -20°C for downstream uses.

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